



CK20



HE staining in scanner view.



Immunostaining, scanner view

Provided sections in this run:

| TISSUE | PROPORTION | INTENSITY |
|-------------------------------|--|-----------|
| Appendix | A strong, distinct cytoplasmic staining reaction of lining epithelial cells more towards the surface in the appendix and focal weak to moderate staining reaction in crypt cells. | Intense |
| Colonic adenocarcinoma | A weak to strong, distinct cytoplasmic staining reaction of the vast majority of neoplastic cells in the colon adenocarcinoma. | Intense |
| Stomach | Vast majority of the foveolar lining epithelial cells are unstained. Occasional focus of intestinal metaplasia shows moderate to intense cytoplasmic staining. This tissue showed non-specific staining in majority of the participant and hence removed from calculations. However, excessive background staining is factored in. The intestinal metaplasia is not represented in all participants and hence that is also discounted. | Weak |
| Urothelial carcinoma | Weak to moderate intensity staining in nearly 60% of the tumor cells. | Moderate |



Performance characteristic of the group:

| | |
|------------------------------|-----------|
| N | 35 |
| MEDIAN OVERALL SCORE | 50 |
| AVERAGE OVERALL SCORE | 40 |

Criteria for assessment:

| | Optimum | Good | Borderline | Poor |
|----------------------|----------------|-------------|-------------------|-------------|
| Overall score | >49 | 40-49 | 36-40 | <36 |

Break up of results

| Result | Participants |
|------------------------------|---------------------|
| Optimum | 26 |
| Good | 3 |
| Borderline | 1 |
| Poor | 5 |
| Section not submitted | 9 |

Break up of results based on technique

| Result Technique | Optimum | Good | Borderline | Poor | Total |
|--|----------------|-------------|-------------------|-------------|--------------|
| Automated | 1 | 0 | 0 | 1 | 2 |
| Automated Ventana Bench Mark XT | 4 | 0 | 0 | 2 | 6 |
| Automated Ventana Benchmark Ultra | 0 | 0 | 0 | 1 | 1 |
| Automated, Ventana Benchmark GX | 1 | 0 | 0 | 1 | 2 |
| Automated, Ventana Medical systems, Inc., Roche Diagnostics | 1 | 0 | 1 | 0 | 2 |
| Manual | 18 | 2 | 0 | 0 | 20 |
| NA | 1 | 1 | 0 | 0 | 2 |



Clones used in the run by participants

| Clone | Participants |
|----------|--------------|
| EP23 | 8 |
| GM003 | 1 |
| IT-KS208 | 2 |
| Ks 20.8 | 6 |
| Ks20.8 | 10 |
| NA | 11 |
| SP33 | 5 |
| ZM42 | 1 |

Break up of results based on the clones used

| Clone | Result | Optimum | Good | Borderline | Poor | Section not submitted | Total |
|----------|--------|---------|------|------------|------|-----------------------|-------|
| EP23 | | 8 | 0 | 0 | 0 | 0 | 8 |
| GM003 | | 1 | 0 | 0 | 0 | 0 | 1 |
| IT-KS208 | | 1 | 0 | 0 | 1 | 0 | 2 |
| Ks 20.8 | | 3 | 1 | 0 | 1 | 1 | 6 |
| Ks20.8 | | 8 | 1 | 0 | 1 | 0 | 10 |
| NA | | 2 | 1 | 0 | 0 | 8 | 11 |
| SP33 | | 2 | 0 | 1 | 2 | 0 | 5 |
| ZM42 | | 1 | 0 | 0 | 0 | 0 | 1 |



Break up of results based on the vendors

| Vendor | Result | Optimum | Good | Borderline | Poor | Total |
|--|--------|---------|------|------------|------|-------|
| Biocare | | 1 | 0 | 0 | 0 | 1 |
| BIOCARE PM062A | | 1 | 0 | 0 | 0 | 1 |
| BiogeneX AM315-5M | | 1 | 0 | 0 | 0 | 1 |
| Biosystems, PDM049 | | 1 | 0 | 0 | 0 | 1 |
| Cell Mourque Cat 320-M-16-RUO | | 1 | 0 | 0 | 0 | 1 |
| Company -Quartett catalogue number- KS20.8 | | 1 | 0 | 0 | 0 | 1 |
| Dako | | 1 | 0 | 0 | 0 | 1 |
| DAKO IS777 | | 3 | 1 | 0 | 1 | 5 |
| DAKO M7019 | | 1 | 0 | 0 | 1 | 2 |
| Master Diagnostica; 005105QD-R-7 | | 1 | 0 | 0 | 0 | 1 |
| NA | | 1 | 1 | 0 | 0 | 2 |
| PATH N SITU , CAT# HAR036 | | 1 | 0 | 0 | 0 | 1 |
| Pathnsitu | | 3 | 0 | 0 | 0 | 3 |
| PathnSitu Catalogue No: MM1153 | | 0 | 1 | 0 | 0 | 1 |
| PathnSitu, PR036 | | 5 | 0 | 0 | 0 | 5 |
| quartett, P-C005-70 | | 1 | 0 | 0 | 0 | 1 |
| Ventana | | 1 | 0 | 0 | 0 | 1 |
| Ventana (Catalogue No: AM3150321) | | 0 | 0 | 0 | 1 | 1 |
| Ventana (240)05587760001 | | 0 | 0 | 0 | 1 | 1 |
| VENTANA 790-4431 | | 1 | 0 | 1 | 1 | 3 |
| Zeta corporation, 180 | | 1 | 0 | 0 | 0 | 1 |



Break up of results based on format

| Result Format | Optimum | Good | Borderline | Poor | Section not submitted | Total |
|------------------|---------|------|------------|------|-----------------------|-------|
| Concentrated | 3 | 0 | 0 | 0 | 0 | 3 |
| NA | 1 | 1 | 0 | 0 | 8 | 10 |
| RTU | 22 | 2 | 1 | 5 | 1 | 31 |

Break up of results based on temperature of dewaxing

| Result Dewax Minutes | Optimum | Good | Borderline | Poor | Section not submitted | Total |
|-------------------------|---------|------|------------|------|-----------------------|-------|
| 37 | 1 | 0 | 0 | 0 | 0 | 1 |
| >75 | 0 | 0 | 0 | 0 | 0 | 0 |
| 55-60 | 7 | 1 | 0 | 1 | 0 | 9 |
| 61-65 | 8 | 1 | 0 | 0 | 0 | 9 |
| 66-70 | 2 | 0 | 0 | 0 | 1 | 3 |
| 71-75 | 4 | 0 | 1 | 3 | 0 | 8 |
| NA | 1 | 1 | 0 | 0 | 8 | 10 |
| RT | 1 | 0 | 0 | 0 | 0 | 1 |



Break up of results based on pH of retrieval buffer

| Result pH | Optimum | Good | Borderline | Poor | Section not submitted | Total |
|--------------|---------|------|------------|------|-----------------------|-------|
| 6.0-6.9 | 2 | 1 | 0 | 0 | 0 | 3 |
| 7.0-8.0 | 3 | 0 | 0 | 1 | 0 | 4 |
| 8.1-9.0 | 17 | 1 | 1 | 4 | 1 | 23 |
| >9.0 | 0 | 0 | 0 | 0 | 0 | 0 |
| NA | 2 | 1 | 0 | 0 | 8 | 3 |


Break up of results based on incubation time of primary antibody

| Result Time in minutes | Optimum | Good | Borderline | Poor | Section not submitted | Total |
|---------------------------|---------|------|------------|------|-----------------------|-------|
| <15 | 0 | 0 | 0 | 1 | 0 | 1 |
| 15-30 | 6 | 0 | 0 | 2 | 0 | 8 |
| 31-45 | 7 | 0 | 1 | 2 | 0 | 10 |
| 46-60 | 12 | 2 | 0 | 0 | 1 | 15 |
| NA | 2 | 1 | 0 | 0 | 8 | 11 |


General observations in Run 17




| VARIABLE | 1 | 2 | 3 |
|------------------------|---|--|---|
| Technique | Automated, Ventana Benchmark GX | Manual Technique | Manual |
| Clone | SP33 | Ks20.8 | GM003 |
| Vendor | VENTANA BENCHMARK GX , H09036NA | DAKO IS777 | PathnSitu, catalogue number: PR036 |
| Format | RTU | RTU | RTU |
| Batch/Year | HAR036,2022 | 20082795 | R08036MA, 2020 |
| Expiry | 2023 | 19/05/2022; validation done and results are satisfactory | 06/23 |
| Dewaxing temperature | 75 degree celcius | 65 degree | 58-60°C |
| Retrieval | HIER | HIER EDTA Buffer pH 9.0 | HIER |
| Enzyme | NA | --- | NA |
| HIER | company system - cell conditioning 1 - CC1 | Micro wave | Pressure cooker |
| Peak T and Time | 95 degrees and 44 mins | 98 degree 20 minutes | NA |
| Peak Pressure and Time | NA | 98 degree 20 minutes | Pressure cooker, 3 whistles, 5 min |
| Retrieval Buffer | Tris EDTA buffer | In house | In house made |
| pH | 7.8 | EDTA Buffer pH 9.0 | pH:9 |
| Blocking | 3% Hydrogen peroxide | 3 % Hydrogen peroxide with methanol 15 minutes | 10-12 minutes |
| Wash sol | TRIS based buffer solution | TBST buffer pH 7.4 | PBS Wash solution |
| Dilution of RTU | NO | --- | NO |
| Dilution of conc | NA | --- | RTU |
| Diluent | NA | 0 | NA |
| Inc time of Primary | 60 mins | 30 minutes | 1 hour |
| Detection | Ultraview universal DAB detection kit | Polymer based system | Poly Excel HRP/DAB Detection system |
| Cat No | Ultraview universal DAB detection kit Ventana | Master Polymer plus detection system (Peroxidase) | Poly Excel HRP/DAB Detection system Pathnsitu |
| Inc time of Sec | 16 mins | 30 minutes | 15 minutes |
| Chro-substrate | 8 mins | 4-6 minutes | 3-5 minutes |
| Post-treatment | copper, 5 mins | --- | No |
| Counterstain | Hematoxylin | Hematoxylin | Hematoxylin |




Focus of intestinal metaplasia in gastric biopsy




The focus stains with CK20. However, many participants failed to show this up, may be due to non-representation of the focus in their sections.




Optimum staining of colonic adenocarcinoma




Suboptimal staining of colonic adenocarcinoma due to inadequate epitope retrieval.




Optimum staining of appendicular crypts. The staining is more pronounced towards surface while the basal crypts also stain with moderate intensity. There is no staining of any other cell.




Suboptimal protocol staining has resulted in huge background staining and non-specific staining of connective tissue cells



Non-specific staining of gastric mucosa due to inappropriate protocol setting.



Optimum staining of urothelial carcinoma.



Optimum staining of urothelial carcinoma